Interactive Automated Bragg Peak Identification with 3D Neutron Scattering Data



Matthew Kilpatrick*, David Bruhwiler, Evan Carlin, Ben Gur, Robert Nagler (RadiaSoft LLC, Boulder, USA);
Christina Hoffmann, Zachary Morgan, Andrei T Savici, Matthew Tucker, Bogdan Vacaliuc (ORNL, Oak Ridge, USA);
Alexander Kuhn, Jörg Mensmann, Peter Messmer, Marc Nienhaus, Steffen Roemer, Dragos Tatulea (NVIDIA, Berlin, Germany)
*kilpatrick@radiasoft.net





Sirepo

Benzil

Cluster ID


100

Summary


- Fast, interactive 3D rendering of reciprocal space
- o in your browser via NVIDIA's IndeX technology
- Automated identification of Bragg peaks
- DBSCAN clustering distinguishes peaks from noise
- o more accurate and robust than Mantid algorithms
- Mantid uses DBSCAN peaks for UB matrix & hkl indices
- o predicts 100% of the indices within 0.15 tolerance
- Mantid-only approach yields 76.5% within tolerance
- Statistical fit on the cluster size can distinguish a complex diffuse background from both noise and Bragg peaks
- Complete analysis of 5 different types of crystal lattices with accurate automated Bragg peak finding


IndeX Visualization

- NVIDIA IndeX technology [3]
 - 3D browser based volumetric interactive framework
 - Scientists can interact with massive data sets
 - Make real-time modifications
 - Navigate to the most pertinent parts of the data in real time
- Our implementation of IndeX
 - Docker based deployment Interactively view many gigabytes of neutron single-crystal diffraction data.



- Live Data Streaming!
 - Successful implementation of IndeX at ORNL with live data
 - Interactive neutron intensity analysis
 - Interactive color mapping







Background

- TOPAZ instrument [I] at the Spallation Neutron Source at Oak Ridge National Laboratory
- Typical analysis using Mantid [4]


SAK RIDGE National Laboratory

- Non-interactive large dataset analysis
- Methods use 2D slices which can miss important features






Bragg Peak Detection with Machine Learning



- Density based clustering algorithm with DBSCAN
- Can identify oblong cluster and attribute data to noise
- Sparse background can be easily removed
- Automatic peak finding within datasets
 Calculate a signal-to-noise for identified peaks
- Tunable cut to allow for optimized data cleaning
- Comparisons with Mantid
 - Compare peak locations with clean signal-to-noise
 - Complete analysis demonstration with 5 different crystal lattices at 10 K and 300 K
 - DBSCAN shows a higher number of peaks within reconstruction tolerance
 - Automatic scanning of signal-to-noise and intensity parameter space to optimize Bragg peak reconstruction
- Further work...
 - Normalization to merge datasets from different user runs
- Development neural network to remove consistent background from detector devices



References

I. L. Coates et al., Rev. Sci. Instrum. 89, 092802 (2018). doi:10.1063/1.5030896

2. K. Bruhwiler et al., IPAC Proceedings, TUPAB413 (2021), https://accelconf.web.cern.ch/ipac2021/papers/tupab413.pdf

3. The IndeX home page, https://developer.nvidia.com/nvidia-index 4. Mantid (2013): Manipulation and Analysis Toolkit for Instrument Data.; Mantid Project. url: http://dx.doi.org/10.5286/SOFTWARE/MANTID.





