## Interactive Automated Bragg Peak Identification with 3D Neutron Scattering Data



Matthew Kilpatrick\*, David Bruhwiler, Evan Carlin, Ben Gur, Robert Nagler (RadiaSoft LLC, Boulder, USA); Christina Hoffmann, Zachary Morgan, Andrei T Savici, Matthew Tucker, Bogdan Vacaliuc (ORNL, Oak Ridge, USA); Alexander Kuhn, Jörg Mensmann, Peter Messmer, Marc Nienhaus, Steffen Roemer, Dragos Tatulea (NVIDIA, Berlin, Germany) \*kilpatrick@radiasoft.net





### Summary

- Fast, interactive 3D rendering of reciprocal space
- In your browser via NVIDIA's IndeX technology
- Automated identification of Bragg peaks
- DBSCAN clustering distinguishes peaks from noise
- More accurate and robust than Mantid algorithms
- Mantid uses DBSCAN peaks for UB matrix & hkl indices
- Predicts 100% of the indices within 0.15 tolerance
- Mantid-only approach yields 76.5% within tolerance
- Statistical fit on the cluster size can distinguish a complex diffuse background from both noise and Bragg peaks
- Complete analysis of 5 different types of crystal lattices with accurate automated Bragg peak finding
- Live Demo of Index visualization
  - Interactive visualization of many crystal lattices
  - Intensity analysis to find Bragg peaks
  - Docker based deployment

# Background • TOPAZ instrument [1] at the Spallation Neutron Source at Oak Ridge National Laboratory Typical analysis using Mantid [4] Non-interactive large dataset analysis Methods use 2D slices which can miss important features

### IndeX Visualization

25 detector installed



#### NVIDIA IndeX technology [3]

- browser based volumetric interactive framework
- Scientists can interact with massive data sets
- Make real-time modifications
- Navigate to the most pertinent parts of the data in real time
- Our implementation of IndeX
  - Docker based deployment Interactively view many gigabytes of neutron single-crystal diffraction data.
- Live Data Streaming!
  - Successful implementation of IndeX at ORNL with live data
  - Interactive neutron intensity analysis
  - Interactive color mapping


#### Live Demo!


- View five different crystal structures on IndeX
- Full 3D capabilities
- Color mapping customization for conference ready visualizations

#### Future Features

- Automatic UB matrix calculation and visualization
- Streaming control and interactability

SAK RIDGE National Laboratory






# **Automated Bragg Peak Detection**

#### Density based clustering algorithm with DBSCAN

- Can identify oblong cluster and attribute data to noise
- Sparse background can be easily removed
- Automatic peak finding within datasets
- Calculate a signal-to-noise for identified peaks
- Tunable cut to allow for optimized data cleaning
- Comparisons with Mantid
  - Compare peak locations with clean signal-to-noise
  - Complete analysis demonstration with 5 different crystal lattices at 10 K and 300 K
  - OBSCAN shows a higher number of peaks within reconstruction tolerance
  - Automatic scanning of signal-to-noise and intensity parameter space to optimize Bragg peak reconstruction
- Further work...
  - Normalization to merge datasets from different user runs
  - Development neural network to remove consistent background from detector devices




#### References

I. L. Coates et al., Rev. Sci. Instrum. 89, 092802 (2018). doi:10.1063/1.5030896

2. K. Bruhwiler et al., IPAC Proceedings, TUPAB413 (2021), https://accelconf.web.cern.ch/ipac2021/papers/tupab413.pdf

3. The IndeX home page, https://developer.nvidia.com/nvidia-index









4. Mantid (2013): Manipulation and Analysis Toolkit for Instrument Data.; Mantid Project. url: http://dx.doi.org/10.5286/SOFTWARE/MANTID.