

Interactive Visualization and Automated Analysis of Neutron Scattering Experiments

Matthew Kilpatrick*, David Bruhwiler, Evan Carlin, Gurhar Kalsa, Boaz Nash (RadiaSoft),

Christina Hoffmann, Zachary Morgan, Andrei T Savici, Matthew Tucker, Bogdan Vacaliuc (ORNL)

Alexander Kuhn, Jörg Mensmann, Peter Messmer, Marc Nienhaus, Steffen Roemer, Dragos Tatulea (NVIDIA)

*kilpatrick@radiasoft.net

March 7, 2024

ICFA 2024, Lahan Select Gyeongju, Republic of Korea



This material is based upon work supported by the U.S. Department of Energy, Office of Science, Office of Basic Energy Sciences DE-SC0021551.



Boulder, Colorado USA | radiasoft.net



NVIDIA




U.S. DEPARTMENT OF
ENERGY

Office of
Science

Neutron Scattering probes within Materials


- Why do we shoot particles at things?
 - X-ray scattering probes low-Z materials and are scattered by electrons
 - Electron scattering is great for low-Z materials and is strongly interacting
- Neutrons are weakly scattering!
 - Neutrons will either pass straight through most materials or only single scatter off the atomic nuclei in the material
 - Neutrons can penetrate materials that stop X-rays or electrons
 - Beneficial for studying materials near absolute zero



<https://www.ornl.gov/blog/what-makes-neutron-scattering-unique>


TOPAZ Single-crystal Diffractometer

- Located at ORNL and receives neutrons from the Spallation Neutron Source
- Samples can be measured with high precision for volumetric sampling in reciprocal space (momentum measurements)
 - Samples are rotated to measure all aspects of the lattice
 - Temperature control from 5 K - 450 K
 - Broad Q coverage






25 detector installed

OAK RIDGE
National Laboratory




<https://neutrons.ornl.gov/topaz>

TOPAZ Single-crystal Diffractometer



Analysis methods for TOPAZ data


- **Mantid**
 - Open source community developed application
 - Algorithms use raw data processed within 2D slices of the 3D data
- **Some challenges faced by the community**
 - Large datasets can be up to 100 Gb in size and current tools limit interactivity
 - Displaying data requires a high level of user interaction
 - 2D slices may miss key features of the diffraction data
 - May be slow to run due to optimization shortfalls
 - Timescale: 1-2 hours
- **Working collaboratively to leverage Mantid**
 - Machine learning – DBSCAN, U-Nets
 - Automated Bragg peak identification
 - Automated background subtraction
- **Also improving the interactivity of the data viz**
 - Close collaboration with the NVIDIA IndeX development team (Berlin)
- **Working towards combination of 3D viz and automated analysis**



<https://www.mantidproject.org/>

Analysis of Scattering Data for a Single Silicon Crystal


- Preprocessing the data
 - Measurements of the particles in momentum space
 - Looking for pockets of constructive interference
 - Measurements projected onto a 3D mesh
- DBSCAN is a density based clustering algorithm
 - Can identify oblong clusters and attribute data to noise if it doesn't match anything
 - Sparse background can be removed easily
 - Provides real time identification of Bragg peaks
 - Agrees with Mantid reconstruction of Bragg peaks



<https://scikit-learn.org/stable/modules/generated/sklearn.cluster.DBSCAN.html#sklearn.cluster.DBSCAN>

Background Removal Processes

- **Sample mounting**
 - Crystal sample can be tested at low temperatures
 - Aluminum cylinders have a randomly oriented crystal structure
 - These rings will scatter neutrons into a sphere or ring
- **Scattered neutrons**
 - Each plane orientation will create rings within reciprocal space
 - Translates into a sphere in 3D representation



TOPAZ Crygoniometer

Sample mount

IVC Shield will thread on
and off, one to two turns

Sample

Sapphire Window
(Optional)

Sapphire
Window


OVC O-ring allows OVC to
slide on and off.
Vacuum will hold it in position

Identify Accurate Ground Truths

- **Properties of Aluminum Rings**
 - Monochromatic in neutron intensity
 - Concentrated at a fixed radius
 - Varies from sample to sample
- **Define a scale factor**
 - Simple clustering
 - Tune DBSCAN parameters to get rings
 - Gradient based definition for clustering
 - Normalized by the mean number of voxels

$$SF_i = \max\left(\sum_i^{N+i} \nabla w_i\right)$$


- How do these two methods compare?




Aluminum Labels

Aluminum 47829

DBSCAN





Aluminum 47829



Identification of Aluminum Rings


- Building a neural network
 - Identify a 3D CNN to train and predict aluminum rings
 - General Unet mainly used for image segmentation
 - Fully connected
 - Contraction and expansion paths
 - 3D Unet used for learning dense volumetric segmentation from sparse annotations






Identification of Aluminum Rings

- Building a neural network
 - Identify a 3D CNN to train and predict aluminum rings
 - General Unet mainly used for image segmentation
 - Fully connected
 - Contraction and expansion paths
 - 3D Unet used for learning dense volumetric segmentation from sparse annotations

Labelled Rings





Full Sample




Identification of Aluminum Rings

- Building a neural network
 - Identify a 3D CNN to train and predict aluminum rings
 - General Unet mainly used for image segmentation
 - Fully connected
 - Contraction and expansion paths
 - 3D Unet used for learning dense volumetric segmentation from sparse annotations


Full Sample




Labelled Rings



Predicted Rings




System Architecture, with JupyterHub Interface



NVIDIA's IndeX: Browser-based 3D visualization

- IndeX: interactive 3D volumetric visualization framework
 - rotate, pan, zoom into large 3D volumes
 - executes remotely on a GPU server
 - streams video to the user's browser
- TOPAZ event data is converted into a voxel database
 - voxels are a regular 3D grid
 - each is colored and/or made transparent based on the number of neutron events within it
 - colors & transparency can be dynamically edited
- The voxel database can be dense or sparse
 - large fraction of voxels are empty, so the sparse data format is very efficient
- Interactive response is robust to WiFi quality
 - even with a poor network connection, the typical delay after a mouse click is < 0.25 s





Summary

- Take the user out of the loop
 - Automated sample identification
 - Real time sample visualization and analysis
 - Information at your fingertips
- 3D visualization
 - Integrated on a single GPU
 - Interactive analysis
- Open-source collaboration
 - Applications to many real time 3D datasets
 - Broad applications to samples within a beamline

THANK YOU!

Acknowledgements

- Our ORNL partners
 - C. Hoffman, B. Vacaliuc, J. Taylor, J. Kohl, Z. Morgan, A. Savici,
- NVIDIA
 - M. Nienhaus, A. Kuhn, D. Tatulea, J. Mensmann, S. Roemer,
- RadiaSoft
 - D. Bruhwiler, E. Carlin, G. Kalsa, B. Nash
- DOE Office of Science, Basic Energy Sciences
 - SBIR and STTR Program, award # DE-SC0021551

Disclaimer

This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.