

# Automation of Alignment and 3D Analysis for Neutron Scattering Experiments



Matthew Kilpatrick  
AI for Scattering – UMD




This material is based upon work supported by the U.S. Department of Energy, Office of Science,  
Office of Basic Energy Sciences, under Award Number DE-SC0021551 and DE-SC0021555.




# TOPAZ Single -crystal Diffractometer

- Located at ORNL and receives neutrons from the Spallation Neutron Source
- Samples can be measured with high precision for volumetric sampling in reciprocal space (momentum measurements)
  - Samples are rotated to measure all aspects of the lattice
  - Temperature control from 5 K- 450 K
  - Broad Q coverage






25 detector installed

 OAK RIDGE  
National Laboratory





# TOPAZ Single -crystal Diffractometer




# Background & Motivation


- Sample alignment
  - Neutron production time is limited
  - Some activities require constant realignment, such as temperature scans
  - User facilities especially face schedule constraints
  - Machine learning (ML) is a key automation tool
- Alignment protocols vary between beamlines
  - Opportunity to employ & test models
  - Broad applications for sample alignment

# UNet Architecture for Image Segmentation


Regular Convolutions  
Expand feature space



Regular Convolutions  
Compress feature space





Concatenations  
Improve spatial precision




Max Pooling

Reduce spatial dimensions



Up-Convolution

Expand spatial dimensions



Feature -space Representation

**Contraction Path** →


Extract features, lose spatial information

**Expansion Path** →

Expand features, regain spatial information


# Machine Learning for Sample Identification

- Supervised learning problem
  - Input images come from optical camera
  - Beamline scientist / expert identifies the sample mask
- Train neural network model to learn the image to mask relationship
- Find the sample CoM






# Uncertainty Quantification



- During supervised training
  - Real error compared to human defined masks
- During testing and operations
  - Ground-truth data not available
  - Employ statistics from ensemble predictions
  - Variance between many trained models
- Takeaways
  - Excellent CoM prediction
  - Negligible ensemble spread with low uncertainty
  - Mask variance is restricted to edges

# Controls Problem

- Pass images through ML model
  - Retrieve predicted mask & center of mass
  - Retrieve uncertainty measurement
- Evaluate quality of predictions/state of controls
- Pass values back to beamline controls
  - Sample/beam center of mass offsets
  - Human intervention needed?
- Retrain networks with feedback from cycle





# Automated Alignment

rocess (press any key)

-----  
-----

Beamline Processes:


lightswitch\_cryo  
align\_cryo  
temp\_ramp

Interface Actions:

- 0) Exit
- 1) Print beamline state
- 2) Print beamline element state
- 3) Time a beamline process


Please choose a beamline processes or interface action:

-> █



# Analysis methods for TOPAZ data


- **Mantid**
  - Open source community developed application
  - Algorithms use raw data processed within 2D slices of the 3D data
- **Some challenges faced by the community**
  - Large datasets can be up to 100 Gb in size and current tools limit interactivity
  - Displaying data requires a high level of user interaction
  - 2D slices may miss key features of the diffraction data
  - May be slow to run due to optimization shortfalls
- **Working collaboratively to leverage Mantid**
  - Machine learning DBSCAN
    - Reduces user interaction
    - Automated Bragg peak identification
- **Also improving the interactivity of the data viz**
  - Close collaboration with the NVIDIA IndeX development team (Berlin)
- **Working towards combination of 3D viz and automated analysis**



<https://www.mantidproject.org/>

# Analysis of Scattering Data for a Single Silicon Crystal

- Preprocessing the data
  - Measurements of the particles in momentum space
    - Looking for pockets of constructive interference
  - Measurements projected onto a 3D mesh
- DBSCAN is a density based clustering algorithm
  - Can identify oblong clusters and attribute data to noise if it doesn't match anything
  - Sparse background can be removed easily



<https://scikit-learn.org/stable/modules/generated/sklearn.cluster.DBSCAN.html#sklearn.cluster.DBSCAN>

# Comparison Cluster Identification with Mantid

- Peak locations

- Fast Fourier Transform (FFT) to calculate UB matrix (orientation matrix) and predict Miller indices

- Limitations of Mantid peak finding

- Needs to be told how many peaks to look for
- Use the number of peaks that is found by DBSCAN

$$\underbrace{Q_s}_{\text{Measured Peaks}} = 2\pi \cdot \underbrace{UB}_{\text{Calculated with FFT from Mantid}} \cdot \underbrace{\begin{pmatrix} h \\ k \\ l \end{pmatrix}}_{\text{Predicted Miller Indices}}$$

- What do we compare?

- UB matrix is not unique
  - Due to the symmetry of the Bragg peaks you can rotate the matrix and still predict the peaks in the proper location
- Miller indices are not unique
  - Due to the rotation of the UB matrix different hkl indices can map to the same peak
  - Confirm values are 'close' to an integer value
- Compare peak locations
  - Clean with Signab-noise

<https://docs.mantidproject.org/v4.0.0/algorithms/FindUBUsingFFT-v1.html>



<https://docs.mantidproject.org/v4.0.0/concepts/Lattice.html#lattice>



# Quantifying the validity of Bragg peaks

- Calculate a Signal-to-Noise for each matched peak
  - A Bragg peak should have a large S/N ratio
  - Background events result in a minimal ratio
  - Tunable cut to allow for optimized data cleaning

$$\frac{S}{N} = \frac{I_{\delta S}}{\sqrt{I_{\delta B}^2 + N_S}}$$






# Silicon 10K and 300K unit cell (reduced) parameters

- Good agreement

- Mantid and DBSCAN methods show the same lattice parameters for Silicon
- Standard deviation is comparable!



- DBSCAN datasets

- Automated peak finding allows for similar results to Mantid only methods
- Files with only DBSCAN results are shifted from the normal results
- All lattice parameters agree with the theoretical value of Silicon




# Silicon 10K and 300K

- Predict UB matrix with identified Bragg peaks
  - DBSCAN finds 100% of Bragg peaks within tolerance
  - Mantid finds 55% of Bragg peaks within tolerance
  - Silicon 10K
    - 24 3D silicon datasets could not find a UB matrix with Mantid peaks
  - Silicon 300K
    - 2 3D silicon datasets could not find a UB matrix with Mantid peaks



# ADARA Live Streaming into NVIDIA's IndeX



# Summary

- Take the user out of the loop
  - Automated sample identification
  - Real time sample visualization and analysis
  - Information at your fingertips
- 3D visualization
  - Integrated on a single GPU
  - Interactive analysis
- Open-source collaboration
  - Applications to many real time 3D datasets
  - Broad applications to samples within a beamline

THANK YOU!

# Acknowledgements

- Our ORNL partners
  - C. Hoffman, S. Calder, B. Vacaliuc, J. Taylor, J. Kohl, Z. Morgan, A. Savici, B. Krishna, R. Gregory & G. Taufer
- NVIDIA
  - M. Nienhaus, A. Kuhn, D. Tatulea, J. Mensmann, S. Roemer,
- RadiaSoft
  - D. Bruhwiler, E. Carlin, J. Edelen, M. Henderson, I. Pogorelov,
- DOE Office of Science, Basic Energy Sciences
  - SBIR and STTR Program, award # DE-SC0021555, DE-SC0021551